


Objectives

- Learn how to make blood smears
- Demonstrate the ability to examine a stained smear under the microscope and accurately perform a manual differential, following the proper procedural steps.
- Identify and differentiate all major blood cell lines, including both mature and immature forms.

Slide Making (Push Smear)




Requirements for Proper Smear Preparation:

- Clean glass slides
- Label frosted side up- label according to SOP
- Small blood drop
- Proper placement of drop
- Smooth spreading of drop
- Peripheral blood smear made from EDTA anticoagulated blood.
- Ideally made within 3 hours


Step 1:

Place a small drop of mixed venous blood on a glass microscope slide.




Step 2:

A spreader slide is positioned at angle and slowly drawn toward the drop of blood.




Step 3:


The spreader slide has been brought in contact with the drop of blood and is being drawn away...



Step 5: Allow slide to dry 5-10 minutes




Types of Stains



Wright's Stain



Leishman Stain



Giemsa Stain



Some stains can look really beautiful, Giemsa for malaria.

Frequent Staining Problems



- Wright's Stain "Too Alkaline"




Too Alkaline Stain:

- 1) thick blood smear
- 2) prolonged staining
- 3) insufficient washing
- 4) alkaline pH of stain components

Correction:

- 1) check pH
- 2) shorten stain time
- 3) prolong buffering time

Frequent Staining Problems



- Wright's Stain "Too Acidic"



Too Acid Stain:

- 1) insufficient staining time
- 2) prolonged buffering or washing
- 3) old stain

Correction:

- 1) lengthen staining time
- 2) check stain and buffer pH
- 3) shorten buffering or wash time


GCLP standards for Staining QC Slides should be:

- Performed daily or day of use
- When a new batch of stain is made or if there is a change in lot number
- Documented on an appropriate QC record form
- Troubleshooting should also be documented as it occurs



In the hand: an example of Slide staining QC scheme. Please note this is only an example and your laboratory must modify and implement a QC scheme base on your laboratory policies and needs.


Examination of Peripheral Blood Smear



<http://www.lab.anhb.uwa.edu.au/mb140/corepages/blood/blood.htm>

What to Look for & Report


- WBC estimates
- Platelet estimates
- Numeration of WBC inclusions, RBC morphology and inclusions
- Manual differential



WBC Estimation

Performance of the leukocyte estimate using the high power (50x) oil immersion objective

$$\text{WBC Estimate (count/mm}^3\text{)} = \text{Average of WBC in 10 field counted} \times 3000$$




Count the number of both intact and disrupted leukocytes in 10 microscopic fields where the erythrocytes are partially but not completely overlapped (zone of morphology). Divide the total number by 10 to establish the mean number of leukocytes per field. Multiply the mean by 3000 to determine the estimated WBC count/mm³. These estimates should approximate that obtained by the cell analyzer, the automated leukocyte count/mm³ is less than 25,000 the estimate number should agree within 20% of the automated count. If the estimate does not agree within 20%, repeat the WBC count on the instrument and/or prepare another smear to ensure that the correct sample was tested and/or used to prepare the smear. This method is less reliable when the automated leukocytes count/ mm³ is greater than 25,000 and your laboratory has to establish policies to verify this high WBC count, these can include repeating the samples on the automated analyzer or performing a dilution.

1. Count the number of WBCs in 10 field
2. Divide the total number by 10 to establish the mean
3. Multiply the mean by 3000 to determine estimated WBC count/mm³.

Platelet Estimation

Performance of the platelet estimate using the high power (100x) oil immersion objective

Platelet Estimate =
Average of platelets in 10 field counted X 20,000



- If a significant numbers of giant platelets and/or platelet clumps are detected, a peripheral estimate of WBCs should be done to prevent reporting spuriously high white blood cell count.

- Do a manual platelet estimate count if a significant numbers of Microcytic red blood cells and/or small cell fragments are detected during the RBC morphology.

1. locate an area of approximately 150 red blood cells
2. Count all the platelets in that area
3. Repeat this until 10 field have been counted
4. Keep the total number of all 10 areas counted
5. Divide the total number by 10 to obtain the average platelet per field, then multiply by (x) 20,000 to obtain the platelet estimate

In the presence of significant platelet clumping this estimate is not accurate. The laboratory has to establish a policy as how to report the platelet count. Suggestions include; redrawing patient's sample in a blue top tube (sodium citrates which is the anticoagulant choice that prevents the in vitro platelet clumping), vortexing the sample (please see you hand out for a copy of this procedure) or only report-platelets appear decrease, adequate or increase and do not provide a platelet count.

RBC Morphology & WBC Review

- ▶ Review a minimum of 10 fields, using the 100x oil immersion objective
 - RBC size and hemoglobin content
 - RBC shapes
 - RBC inclusions
 - WBC inclusions and Toxic Granulation

1+	2+	3+	4+
Occasional/Few	Moderate	Many/Numerous	For toxic granulation
Slight	Moderate	Marked	




The RBC Morphology and WBC Review will be performed with every manual differential or as a separate examination. The test will be performed using the high power (100X) oil immersion objective. Begin the examination in an area where the red cells are touching but not overlapping. There are approximately 300 RBC's present in this area in patients with a normal hematocrit. Locate a mature (normal and not atypical) lymphocyte for comparative purposes and determination of Normocytic, Microcytic or macrocytic RBC.

Each morphologic abnormality observed should be quantitated ("graded") separately as to severity ("slight to marked" or "1+ to 3+"). 4+ is usually for toxic granulation. SOP for grading.

There is a SOP on your flash drive that give you an example how to set up for grading the different features

Differential

- Perform a 100 cell differential count using 50x or 100x objective
- 100x objective must be utilized to classify any immature cells



Picture description: Scanning technique for peripheral blood differential count and morphologic evaluation. (a) Ten microscopic fields are examined in a vertical direction from bottom to top (or top to bottom). (b) The slide is horizontally moved to the next field (c) Ten microscopic fields are counted vertically. (d) The procedure is repeated until 100 leukocytes have been counted (for a 100-cell count). Or Use a differential cell counter, count and classify 100 white blood cells using the cross-sectional technique where the white cells are counted in consecutive fields as the blood film is moved from side to side as pictured

A differential count of at least 100 white blood cells is performed. There are situations when you want to use more or less cells based on your laboratory policy and also the WBC count.



If the automated WBC count is less than $1.0 \times 10^3/\mu\text{L}$, a 50 cell count can be performed and converted to 100% by multiplying all values by 2. Add the following comment in result slip: 50 cells counted during the manual differential and converted to percentage.

If the automated WBC count is greater than $30.0 \times 10^3/\mu\text{L}$, a 200 cell count must be performed and converted to 100% by dividing all values by 2. Round accordingly and report only whole numbers. Add the following comment in the result slip: 200 cells counted during the manual differential and converted to percentage.

Differential- Correction

- Presence of NRBC or Megakaryocytes


$$\frac{\text{WBC} \times 100}{\# \text{ NRBCs and/or Megakaryocytes} + 100} = \text{Corrected WBC}$$




If nucleated red blood cells (NRBC) or Megakaryocytes nucleus are seen during the differential count, enumerate them separately from the white blood cell count.

Your laboratory has to establish a policy as when to use this formula for example: the WBC count requires correction if more than 10 NRBC or Megakaryocytes are counted during the differential count


Blood Cell Overview



<https://www.cdc.gov/parasites/blood.html>



<https://www.quora.com/What-are-the-main-functions-of-white-blood-cells>



<http://eclinpath.com/hematology/morphologic-features/platelets/normal-platelets/bovine-platelets/>

RBC Morphology Categories

- Hemoglobin content
- Size
- Shape
- Inclusions



When performing RBC examination on the blood slide there are 4 categories you should evaluate and they are


Hemoglobin content

Size

Shape

Inclusions.

We will look at the Hgb content first



For hemoglobin content there are 3 possibilities and the way you differentiate them is the amount of central pallor in them.

If no more than 1/3 the diameter of the cell it is called a Normochromic RBC

If more than 1/3 the diameter of the cell it is called a Hypochromic RBC. Usually seen in iron deficiency anemia.

And finally in middle is a polychromic RBC which is slightly bluer in color and usually has no middle pallor

Can anyone tell me what test is performed to evaluate the % of polychromic RBC in the bloodstream?

Reticulocyte count


Poikilocytosis

- This term is used to describe variation in the RBC Shape.
- Never call poikilocytosis unless you also note an abnormal shape.



The next evaluation is the variation of shape of the RBCs which is called Poikilocytosis. The most important thing to remember is if you call poikilocytosis **you must** also include the abnormal shape. You should never report just poikilocytosis. If you do it is telling the doctor "I see something and I am not telling."

There are several shapes that can occur; we will now go over a few of the common ones seen.




On the left are Stomatocytes that look very similar in appearance to a normal RBC the only difference in appearance is that the central pallor is more slit-like not round. Most commonly seen in hereditary stomatocytosis which is a membrane defect. This is also a autosomal dominant inheritance trait making up 10-50% of the RBC population. They are also seen if the blood has a high acidic pH. In the osmotic fragility test these cells will lysed quicker.

Then on other side we have the ovalocytes/elliptocytes. They usually have a rod shape with concentration of hemoglobin on the ends. In normal blood they make up less than 1% of the population. They can be seen in many diseases the most common is Heredity Elliptocytosis which is caused by an abnormality of the RBC skeletal membrane proteins. In this condition they will make up more than 25% of the cell population.

There is speculation that it may be beneficial to have this disorder in malaria prone locations because it seems to have a resistant to it. They can also be seen in sickle cell anemia, severe iron deficiency anemia, thalassemias

What do we have?



So let's see how many shapes you can find on this slide –

Elliptocytes/ovalocytes

Teardrops

Spherocytes

Schistocytes

Macrocytes

Hypochromic

Anisocytosis – variability of the cells in size

Poikocytosis – variability of shape

So you see if you just called Poikocytosis and nothing else the doctor would not know anything.

This patient has Hereditary Pyropoikilocytosis which is a subtype of Hereditary Elliptocytosis -Congenital hemolytic anemia –It is both a defect and deficiency in the spectrin in the cell membrane.


The prefix Pyro in greek means fire

It is named this because in the osmotic fragility test these cells will lyse at the lower temperature than normal


You can also see a similar morphology in a patient that had severe thermal burns.

So not let's look at some inclusion that can be seen with RBC's

Cold Agglutination



Rouleaux



And I just wanted to show you these 2 slides that also give you some problems on Cell ID. Agglutination and rouleaux .

Rouleaux. Is described as stacks of coins. If you see this only in thick area of the slide it is an artifact If it is seen in the thin area then it is significant to note on your report. Most common disease state is multiple myeloma.

Agglutination is the clumping of RBC seen with cold agglutinins most commonly caused by an IgM antibody. It can occur after a viral or Mycoplasma infection, chronic idiopathic cases and underlying lymphoproliferative disorders and PCH paroxysmal cold hemoglobinuria. For the agglutination that occur after a viral infection it is usually limiting that eventually goes away in a couple of months..


Breakout for unknown RBC grading.

RBC Inclusions

Pertains to any intracellular material
found in the RBC



Inclusions pertains to any intracellular material found in the cell




On the left is Basophilic stippling - Basophilic is due to staining aggregates of ribosomes containing RNA and may be seen in thalassemia, lead poisoning and sideroblastic anemias and sickle cell anemia. When it is coarse stippling it is never normal and should be noted.

On the side are Pappenheimer bodies which are composed of iron that are visible with either Wright stain or Iron stain such as Prussian Blue. They can be seen in a variety of disease states including sideroblastic anemia and following a splenectomy, thalassemias, and megablastic anemias.

White Blood Cells


- There are five morphologic types of WBCs normally seen in the peripheral blood.

- Neutrophils
- Lymphocytes
- Monocytes
- Eosinophils
- Basophils



There are five different types of WBC seen normally in the blood and they are Neutrophil, Lymphocytes, Monocytes, Eosinophils and Basophils.

We will look at the granulocytes first



Here we have the maturation stages of the Neutrophil

First four stages are normally seen in the bone marrow. Stage one is the Myeloblast. It is a large cell with a large and round nucleus with delicate chromatin and nucleoli present. You have scant blue cytoplasm around the nucleus with no granules.

Stage 2 is the Promyelocyte – It is similar to the blast but the nucleus is starting to get slightly coarse and the primary granules appear in the cytoplasm.


The third stage is the Myelocyte. The cell is decreasing in size and the nucleus in getting coarser without any nucleoli present. The cytoplasm is changing color it is more light blue to pink. At this stage the secondary granules appear depending on the cell line. For neutrophils they are fine blue pink in color.

Stage 4 is the Metamyelocyte. The nucleus is getting indented or kidney shaped. The chromatin of the nucleus is coarser and the cytoplasm is a light pink with fine pink-purple granules.

Stage 5 is the Banded neutrophil. The nucleus in indented more and the chromatin is getting coarser and darker. They cytoplasm is now light pink with fine pink-purple granules. These can be seen in low numbers in the peripheral blood

The most mature is the segmented neutrophil. It usually contains 2-5 lobes joined by thin filaments with very coarse chromatin. The cytoplasm is pale pink with fine pink-purple granules.

Auer Rods




On the left is a Myeloblast with an Auer rod that can be seen singly or in clusters. The can be seen in myeloid leukemia blast. It is generally considered most diagnostic clue for acute non-lymphocytic leukemia- AML or AMML when you see these.

On the side on top is a NRBC . It has the condensed nucleus and pink cytoplasm.

Below we have a Necrotic or degenerated neutrophil also called a pyknotic cell. This is a neutrophil that has degenerated. The nucleus become dense and homogeneous it can have several lobes or be just one lobe as seen here which could easily be mistaken for a NRBC.


It is important to differentiate these two cells. The nucleated RBC cytoplasm will be pink in color with no granules while Pyknotic neutrophil will have granules.

Pyknotic Cell vs NRBC





Pyknotic cell chromatin condensation. Degenerative condition of a cell nucleus marked by clumping of the chromosomes, hyperchromatism and shrinking of the nucleus

Hypersegmented



Hyposegmented



On the left we see have Hypersegmented Neutrophil that can be seen in B12 and folate deficiencies. This term is used when the cell contains 6 or more lobes. In these anemias they will be seen in greater than 5%.

On the right we have Hypo-segmented neutrophil

The upper slide is of a typical cell seen in Pelger-Huet Anomaly which is an inherited disorder. It usually has the classic bi-lobed nucleus connected by a fine filament making it look like a eyeglass or pince-nez


Beneath is we have the Pseudo-Pelgroid cells.

These cells are form due to abnormal maturation of both the nucleus and cytoplasm. It can cause the nucleus not to mature and divide into lobes but stay a single dense nucleus that could be mistaken for a myelocyte.

VERY IMPORTANT to differentiate that these are not myelocytes but pseudo pelgroid cells. If you call them myelocyte you could lead the doctor to think that they have an immature picture - possibly leukemia or a very bad bacterial infection. They should be given their own category and specified as Pseudo-Pelgroid cells.

These can cells be seen in myelodysplastic syndrome and patients receiving chemo.


Hypogranulation



Here we have Hypogranulation with the lack of granules.


This is usually seen in myeloproliferative diseases but you can also see it in patients with HIV.

Lymphocytes



Lymphoblast

Prolymphocyte




Lymphocyte



For the lymphocyte there was three maturations stages.


First we have the lymphoblast. There are actually 3 different types of lymphoblast. To differentiate them you need to do special testing such as flow to look for cell markers. They usually are a large cell with fine chromatin in the nucleus and several nucleoli. Like all blast there is scant amount of dark blue cytoplasm. They are usually only seen in the bone marrow but if in the blood most likely seen in young children with ALL. Next is the Prolymphocyte. It is a little smaller, condenser nucleus with only one prominent nucleoli. Again not usually seen in peripheral blood but if seen most likely an older adult with CLL. They usually make up about 10% of the population. Then the last stage is the mature lymphocyte. They can be seen both in the Bone marrow and the peripheral blood. They are small compact cells with a round or oval nucleus with very condensed chromatin and no nucleoli. There is small amount of dark blue cytoplasm.

Large Granular Lymphocyte




Some lymphocytes that have more abundant cytoplasm as seen here can contain a few, coarse, unevenly distributed, granules. These are referred to as **large granular lymphocytes**. They are seen in normal blood smears but may be increased in patients with reactive lymphocytes.

Cleft Lymphs




Sézary Cells




On the left we have cleft cells or as my husband calls them the “butt-cracked cells”. Depending on the cause these can be seen in certain infections such as whooping cough others can be seen certain types of lymphomas.

On the right we have Sezary cells. These are found in patients with mycosis fungoides which is a cutaneous T cell lymphoma. They are usually round to oval with folded, grooved, or convoluted nuclear membranes that give them a cerebral-form appearance. The cytoplasm is usually pale, blue to gray and can have many small vacuoles adjacent to the nucleus giving it a “pearl necklace appearance.”



On the left we have Hairy-Cell Leukemia lymphs. These are seen in a chronic lymphoproliferative disease of the B-Cell origin called Hairy Cell Leukemia. It is normally seen in middle-aged elderly males. The classic hairy cell is round to ovoid lymphoid cell that are slightly larger than normal mature lymphs. The usually contain more cytoplasm that is pale blue to grayish blue and the border have the characteristic elongated, fine-hairy cytoplasmic projections. To me it looks like they are trying to push the other cells away. The nucleus is usually finer and can have folds in it. What is the special stain you use to differentiate these cells – TRAP - Tartrate Resistant Acid Phosphatase Stain

On the left frame we have a Basket cell or smudge cell. This name is used for cell death most commonly associated with cells that are fragile and easily damaged in the process of making a blood smear. The cell is usually a lymphocyte. You end up with the nucleus that is either smudged or the chromatin strands spread out from a condensed nuclear remnant giving it the basket-like appearance. Usually there is no cytoplasm or very indistinct. They are associated in some disease states characterized by lymphocyte fragility such as infectious mononucleosis and CLL. One way to avoid these on smears is adding a drop of 22% bovine serum albumin to four or five drops of blood before making the slide. The albumin will cushions the cells and prevent the fragments. Sometimes you can also get fragmented cells on pediatric slides due to trying to spread the thicker blood they may have due to a high Hemoglobin. You end up putting too much pressure on the cells and cause them to fragment. So try using the albumin technique to get a better slide
So let's go onto the Monocyte line




On the left we have large or giant platelets. For Cell ID proficiency that have a normal MCV you should call it a large/giant platelet when the platelet is larger than the average RBC. They are seen in myeloproliferative disorders and myelodysplastic syndromes.



On the other side are hypogranular or agranular platelets. As its name implies it has reduced or absent granularity. They may be normal in size or enlarged and/or misshapen. The cytoplasm stains pale blue or blue gray. They can be associated with myeloproliferative disorders and acute leukemia's or CML

QUIZ TIME

Can you Identify the Cell?





OWA SLIDES



Cell on left- Neutrophil seg/band with Dohle and toxic granulation
Cell on right- lymphocyte



OWA Slides



Cell on left-platelet



Cell on right- monocyte with pseudopod

OWA Slides



Cell on left- Howell Jolly Body
Cell on right- Basophil


OWA Slides



Cell on left- neutrophil seg/band with Dohle and toxic granulations
Cell on right- Band

Blood Cell Identification – Graded

BCP-14



Name this cell: size 15um


- 1. Monocyte
- 2. Abnormal/Reactive Lymph
- 3. Plasma Cell
- 4. Blast

COLLEGE of AMERICAN
PATHOLOGISTS
Laboratory Quality Solutions



Blood Cell Identification – Ungraded

BCP-09



Name this cell: size 20um


1. Monocyte
2. NRBC
3. Blast
4. Abnormal/Reactive Lymph

COLLEGE of AMERICAN
PATHOLOGISTS
Laboratory Quality Solutions



Blood Cell Identification – Ungraded

BCP-09



COLLEGE of AMERICAN
PATHOLOGISTS
Laboratory Quality Solutions

Why is this a Blast:

1. Round to oval 10-20 um
2. N:C Ratio varies- 7:1 - 5:1
3. Large nuclei that has fine lacey chromatin and at least one prominent nucleoli
4. Typically agranulor cytoplasm




Huge nucleus

At least one prominent nucleoli

Very difficult to classify a blast- generally if auer rods are present it is in the myeloid lineage (neutrophil)

Blood Cell Identification – Ungraded

BCP-27



Name this cell: 20um


1. Band
2. Lymphocyte
3. Monocyte
4. Segmented Neutrophil

COLLEGE of AMERICAN
PATHOLOGISTS
Laboratory Quality Solutions



Blood Cell Identification – Ungraded

BCP-27




COLLEGE of AMERICAN
PATHOLOGISTS
Laboratory Quality Solutions

Why is this a Monocyte:

- 1. Round 12-25 um (sometimes have pseudopod extensions)**
- 2. N:C Ratio varies- 4:1 - 2:1**
- 3. Band-like nuclei that has condensed chromatin and absent nucleoli.**
- 4. Blueish gray cytoplasm ground glass appearance**



12-20 Micrometers




oneworld
ACCURACY™

Name this cell:

1. Metamyelocyte
2. Band/Segmented neutrophil
3. Lymphocyte
4. Monocyte


pSMILE



oneworld
ACCURACY™

Name this cell:

- 1. Lymphocyte**
- 2. Band/Segmented Neutrophil**
- 3. Monocyte**
- 4. Blast**



Why is this a Monocyte :

- 1. Round 15-25 um**
- 2. N:C Ratio varies- 4:1 to 2:1**
- 3. Nuclei that has condensed chromatin and absent nucleoli.**
- 4. Blueish gray cytoplasm ground glass appearance with vacuoles**

Acknowledgements

The presenter would like to thank:


- Johns Hopkins University
- Dr. Lori Sokoll- Principal Investigator
- NIH/DAIDS- Daniella Livnat
- Mark Swartz- Project Manager

This project has been funded in whole or in part with Federal funds from the Division of AIDS (DAIDS), National Institute of Allergy and Infectious Diseases, National Institutes of Health, Department of Health and Human Services.

Contract Number: 75N93020C00001



Questions



References

- Reference Leukocyte Differential Count, NCCLS Document H20-A, Vol. 12 No. 1, March 1992.
- Brown, Barbara A., Hematology: Principles and Procedures, Lea and Febiger Book Publisher, Sixth Edition, 1993, Pages 102 to 105.
- Lee, Richard G., Wintrobe's Clinical Hematology; Lea and Febiger Book Publisher, Ninth Edition, 1993, Pages 223 to 238.
- Hoffman, Ronald, Hematology: Basic Principles and Practice, Second Edition, Churchill Livingstone Inc., 1995, Pages 308 to 312.
- Waters, Jerry, Standardization of Red Cell Morphology Reporting Video, CLE (Clinical Laboratory Education), Milwaukee, WI.
- Diggs, L.W., The Morphology of Human Blood Cells Color Atlas, Abbott Laboratories, Inc., Fifth Edition, 1985.
- Glassy, Eric F., Color Atlas of Hematology, College of American Pathologists, 1998.
- <http://www.pathology.vcu.edu/education/PathLab/pages/hematopath/pbs.html>.