



- Two levels of quality control must be run each day that testing is performed.
- It is recommended that testing should be performed by at least 2 different testing personnel.

d. Acceptability criteria:

- pH should be within +/- 1 of the expected value.
- Specific gravity should be within +/- 0.005 of the expected value.
- Semi-quantitative tests should be evaluated based on the comparison method used.
 - if using quality control material, refer to the manufacturer’s acceptable range of each level of quality control.
 - If using EQA, refer to the provider’s participant summary report for acceptability limits.
- For each qualitative or semi-qualitative analyte, use the contingency table below that compares the results of a qualitative test with the outcome of the diagnostic accuracy criteria. The entry in each cell of the table represents the number of specimens corresponding to the labels in the margins.

Method being Validated	Diagnostic Sensitivity and Specificity (Results from Comparison Study)		Total
	Positive	Negative	
Positive	# true positive (TP)	# false positive (FP)	TP+FP
Negative	# false negative (FN)	# true negative (TN)	FN+TN
Total	TP+FN	FP+TN	N

- Calculate the estimated Diagnostic Sensitivity
(True positive rate) = $100 \times [TP/(TP+FN)]$
- Calculate the estimated Diagnostic Specificity
(True negative rate) = $100 \times [TN/(FP+TN)]$
- Calculate the percent Positive Agreement
(Positive Predictive Value)= $100 \times TP/(TP+FP)$
- Calculate the percent Negative Agreement
(Negative Predictive Value) = $100 \times TN/(TN+FN)$
- Compare the results calculated above with the manufacturer’s stated claims for Sensitivity, Specificity and Agreement found in the test kit package insert.

