FINAL VALIDATION REPORT - *M. tuberculosis* culture in LJ media slants

1. Background

Tuberculosis (TB) is the cause of death for 1.5 million people a year [1] and now possibly, we have an increase of affected people [2]. The National TB Control Programs carry out various actions for diagnostic and therapeutic approaches [3] and the TB treatment is monitored through sputum smear microscopy and periodic cultures [3, 4] with the objective of know the exactly category of treatment outcome [5]. While smear microscopy is the most popular method to detect mycobacteria in clinical specimen, culture of the etiological agent remains the accepted “gold standard” for diagnosing mycobacterial infections [6]. Cultivation on solid media, such as LJ, has a critical role in labs. This method is very frequent because has been designed for the isolation and growth of *M. tuberculosis*. Additionally, it is used for the determination of antibiotic susceptibility and microbiological research. [7]. Given some limitations, commercial products are preferred for improving performance; thus, minimizing variability between assays [8].

Since SES participates in clinical trials and this method is useful for diagnostic and treatment monitoring of patients; this report documents the validation process of the BD LJ media slants in our TB lab by comparison with the previously validated LJ media.

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4. Methodology

4.1. Specimens

Spiked simulated samples (sterile water) with standard strains were used for the precision assessments.

- As negative control was used *S. aureus* (ATCC 25923) The spiked specimen was performed using an individually independent bacillar suspension of this strains at McFarland 0.5 with a dilution to 10^{-4}
- As positive control was used *M. tuberculosis* H37Ra (ATCC 25177). The spiked specimen was performed using an individually independent bacillar suspension of this strains at McFarland 0.5 with a dilution to 10^{-5}

- Additionally, we used a blank or system control using phosphate buffer 0.068 M; pH 6.8

Raw available samples used in this validation were chosen from our cryobank (-80°C) with the following characteristics:

- 15 or more specimens with negative culture result by the standard method.
- 15 or more specimens with positive culture result by the standard method.

4.2. Materials and methods

The *M. tuberculosis* culture procedure has performed using:

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4.3. Precision

Precision is the agreement of the measurements of replicate runs of the same sample. It is the process of determining the range of random error.

- 4.3.1. For short-term, *M. tuberculosis* H37Ra, *S. aureus* as non-mycobacteria species and a negative control containing only phosphate buffer were used in triplicate run (individually independent bacterial suspensions), by 2 operators in the Lab.
- 4.3.2. For long-term, *M. tuberculosis* H37Ra, *S. aureus* as non-mycobacteria species and a negative control containing only phosphate buffer were used in triplicate run (individually independent bacterial suspensions), by 2 operators and for three days in the Lab.

Acceptability criteria: Results are approved if the determined reproducibility is 100%.

4.4. Accuracy

Accuracy is the true value of a substance being measured and this verification is the process of determining that the test system is producing true, valid results.

- 4.4.1. At least 30 raw samples were used. These samples were chosen by the results on ZN smear and Xpert MTB/RIF Ultra.

Acceptability criteria: Results are approved if the determined grown is obtained between the 90 - 100% of the specimens.

4.5. Time to Detection (TTD)

For precision and accuracy, TTD was registered for each specimen. This assessment is not applicable, because the specimens used were frozen and the characteristics will be compromised.

4.6. Analytical Sensitivity

This assessment is not applicable, because the specimens used were frozen and the characteristics will be compromised.

4.7. Method Approval

The final decision on methodology validation and acceptance is made by the Laboratory Director after review of all the studies performed as part of the complete method validation process.

4.8. Operators

OPE 01: MUF - Mariella Uchima (Biologist)
 OPE 02: DAO – Diana Alonso (Biologist)

4.9. Supervisor

Nadia Barreda – Biologist, MSc.

5. Results

5.1. Short term (with-in run) precision

5.1.1. Results

Table N° 1. Short term precision results of culture performed on BD LJ Medium Slants using spiked specimens.

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* Bacterial suspension individually prepared for this test. Lot #: Lot number of the LJ slant used for this validation. NG: No grow. COL: Colonies observed, not conclusive. NA: Not applicable.

5.1.2. Assessment of short term (with-in run) precision

Operator MUF 100 %
 Operator DAO 100 %

5.1.3. Acceptability of short term (with-in run) precision

Operator MUF APPROVED
 Operator DAO APPROVED

5.2. Long term (between-day) precision

5.2.1. Results

Table N° 1. Long term precision results of culture performed on BD LJ Medium Slants using spiked specimens.

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* Bacterial suspension individually prepared for this test. Lot #: Lot number of the LJ slant used for this validation. NG: No grow. COL: Colonies observed, not conclusive. NA: Not applicable.

5.2.2. Assessment of long term (between-day) precision

Operator MUF 100 %
 Operator DAO 100 %

5.2.3. Acceptability of long term (between-day) precision

Operator MUF APPROVED
 Operator DAO APPROVED

5.3. Accuracy

Table N° 2. Accuracy results of culture performed on BD LJ Medium Slants using frozen available specimens.

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* Bacterial suspension individually prepared for this test. Lot #: Lot number of the LJ slant used for this validation. NG: No grow. COL: Colonies observed, not conclusive. NA: Not applicable.

5.3.1. Assessment of accuracy precision

Operator MUF 100 %
 Operator DAO 100 %

5.3.2. Acceptability of accuracy precision

Operator MUF	APPROVED
Operator DAO	APPROVED

6. Final considerations

Expectations for culture outcomes, although not perfect, may be reasonable. It is important to note that the samples were chosen for their initial results of ZN smear and Xpert MTB/RIF Ultra, which were performed relatively close to the collection and with fresh samples. Our culture results were obtained for processing of frozen specimens and the bacilli distribution possibly are not the same as fresh samples.

Likewise, it must be considered that although the samples analyzed were raw sputum, they were frozen and it is possible that the thawing process has hydrolyzed the sample and especially the bacilli, with the possibility of viability loss during this process.

However, the results obtained have been quite consistent. We have not found significant differences in the capacity to recover bacilli in the samples analyzed and we hope that the change will not generate any change in the performance of the method in the diagnosis of TB.


With all this, a great ability to practically detect the available positives and an almost perfect similarity between the repetitions carried out by each analyst has been demonstrated, so we consider that the findings are satisfactory to establish this procedure as acceptable.

7. Method Approval

For all the above, the auramine staining method developed by both operators in the SES Lab shows values of reproducibility, accuracy and acceptable detection limits.

FINAL RESULT: *M. tuberculosis* culture method using BD BBL Lowenstein-Jensen Medium Slants REF: 220909 APPROVED.

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